

Home Oxygen Equipment And Consumables 6lpm And Above



baywater
H E A L T H C A R E



High Flow Concentrator

Nidek and AirSep
Intensity - up to 8lpm



Freedom[®] 400 Cylinder **Durations (approx)**

6lpm: 1hour 11mins
(with conserver
3hrs)

8lpm: 53 mins
(with conserver
2hrs 20mins)

10lpm: 43mins

15lpm: 28mins



B10 Cylinder

Durations (approx)

6lpm: 5hrs 53mins

8lpm: 4hrs 25min

15lpm: 2hrs 21min



Conserver **Flow Range**

High Flow Conserver

1 - 8lpm

Standard Conserver

0.5 - 7.5lpm



Cluster Headache Demand Valve

Delivers oxygen as patient breathes, high pressured oxygen 15lpm prescription required



High Flow Stroller Liquid Oxygen Flask Durations (continual)


6lpm: 2hrs 40mins

8lpm: 2hrs

10lpm: 1hr 40mins

12lpm: 1hr 20mins

15lpm: 50mins



Consumables

Cannula Oxygen Concentration

35%: 6lpm

40%: 8lpm

1600Q Cannula: 5-8lpm

1600HF Cannula: 8-15lpm



Consumables

Venturi Mask Oxygen Concentration

31%: 6lpm **Orange**

35%: 8lpm **Yellow**

40%: 10lpm **Red**

60%: 15lpm **Green**



Consumables

Medium concentration masks 2 - 8lpm, approx oxygen concentration

6lpm: 40%

8lpm: 50%

10lpm: 60%

(Advised to assess patients on this mask)



Consumables

Non-rebreathe mask

8 - 15lpm

(approx 80% oxygen concentration)



For More Information

Clinician Hotline:

01270 218050

Portal:

<https://goo.gl/XscMoN>

Website:

www.baywater.co.uk



baywater
H E A L T H C A R E